(EN) INSTRUCTION MANUAL RESCUE STRETCHER I IMMOBILIZER VEST - According to the Regulation (EU) 2016/425	5	(IT) MANUALE DI ISTRUZIONI BARELLA DI SALVATAGGIO I GILET IMMOBILIZZATORE - Conforme al Regolamento (UE) 2016/425	17-18
(CS) NÁVOD K POUŽITÍ ZÁCHRANNÁ NOSÍTKA I IMOBILIZÉROVÁ VESTA - Dle nařízení (EU) 2016/425	5-6	(DE) BEDIENUNGSANLEITUNG RETTUNGSTRAGE I WESTE MIT WEGFAHRSPERRE - Gemäß der Verordnung (EU) 2016/425	18-19
(HU) HASZNÁLATI UTASÍTÁS MENTŐHORDOZÓ I INDÍTÁSGÁTLÓ MELLÉNY - Az (EU) 2016/425 rendelet szerint	6-7	(NL) HANDLEIDING REDDINGSTRETCHER I STARTBLOKKERING VEST - Volgens Verordening (EU) 2016/425	19-20
(PL) INSTRUKCJA OBSŁUGI NOŻE RATUNKOWE I KAMIZELKA IMMOBILIZERA - Zgodnie z Rozporządzeniem (UE) 2016/425	7-8	(ET) KASUTUSJUHEND PÄÄSTEKANDAM I IMMOBILISERI VEST – Määruse (EU) 2016/425 kohaselt	20-21
(RO) MANUAL DE INSTRUCȚIUNI TARGA DE SALVARE I VESTĂ DE IMOBILIZARE - Conform Regulamentului (UE) 2016/425	8-9	(LV) INSTRUKCIJAS Glābšanas nestuves I Imobilizera VESTE - Saskaņā ar Regulu (EU) 2016/425	21-22
(SK) NÁVOD NA POUŽITIE ZÁCHRANNÉ NOSIČE I IMOBILIZÉROVÁ VESTA - Podľa nariadenia (EÚ) 2016/425	9-10	(LT) INSTRUKCIJŲ VADOVAS GELBĖJAS NEŠTUKAS I Imobilizierio liemenė – Pagal reglamentą (EU) 2016/425	22-23
(DA) INSTRUKTIONSMANUAL REDNINGSBÅRE I IMMOBILISER VEST - I henhold til forordning (EU) 2016/425	10-11	(RU) ИНСТРУКЦИЯ ПО ЭКСПЛУАТАЦИИ СПАСАТЕЛЬНЫЕ НОСИЛКИ И ИММОБИЛАЙЗЕРНЫЙ ЖИЛЕТ - В соответствии с Регламентом (EU) 2016/425	23-24
(NO) BRUKSANVISNINGEN REDNINGSBÅRE I IMMOBILISER VEST - I henhold til forordning (EU) 2016/425	11-12	(BG) ИНСТРУКЦИЯ ЗА УПОТРЕБА СПАСИТЕЛНА НОСИЛКА I ЖИЛЕТКА С ИМОБИЛАЙЗЕР - Съгласно Регламента (EU) 2016/425	24-25
(SV) BRUKSANVISNING RÄDDNINGSBAR I IMMOBILISERINGSVÄST - Enligt förordning (EU) 2016/425	12-13	(EL) ΕΓΧΕΙΡΙΔΙΟ ΟΔΗΓΙΩΝ ΔΙΑΣΩΣΤΙΚΟ ΦΟΡΕΙΟ I ΓΙΛΕΚΟ IMMOBILIZER - Σύμφωνα με τον Κανονισμό (EU) 2016/425	25-26
(FI) KÄYTTÖOPAS PELASTUSPAARI I Ajonestottimen liivi - Asetuksen (EU) 2016/425 mukaisesti	13-14	(HR) UPUTSTVO ZA UPOTREBU NOSILKA ZA SPAŠAVANJE I IMOBILIZATOR PRSLUK - Prema Uredbi (EU) 2016/425	26-27
(PT) MANUAL DE INSTRUÇÕES MACA DE RESGATE I COLETE IMOBILIZADOR - De acordo com o Regulamento (UE) 2016/425	14-15	(MT) MANWAL TA' ISTRUZZJONI STRECCH TA' SALVAT I IMMOBILIZER VEST - Skont ir-Regolament (UE) 2016/425	27-28
(ES) MANUAL DE INSTRUCCIONES CAMILLA DE RESCATE I CHALECO INMOVILIZADOR - Según el Reglamento (UE) 2016/425	15-16	(SL) NAVODILA REŠEVALNA NOSILA I JOVNIK ZA IMOBILIZACIJO - Po Uredbi (EU) 2016/425	28-29
(FR) MANUEL D'INSTRUCTIONS BRANCARD DE SAUVETAGE I GILET ANTIDÉMARRAGE - Conformément au Règlement (UE) 2016/425	16-17	(AR) تعليمات يدوي نقالة الإنقاذ سترة منع الحركة - وفقاً ل اللائحة (الاتحاد الأوروبي) 524/6102	29-30


READY FOR:

FALL SAFE

INSPECTOR®

50-7091_STD2023_PO


LISTED ORGANIZATION FOR: (EU) 2016/425 ISO9001:2015

SGS FIMKO OY, P.O. Box 30 (Särkiniementie 3) 00211 Helsinki - Finland
 NOTIFIED BODY: CE 0598

NOTIFIED BODY - CE1019

LISTED ORGANISATION FOR EU TYPE APPROVAL:

VVUÚ. a.s, Pikartská 1337/7 - 716 07 Ostrava - Radvanice, Czech Republic

MEDICAL DEVICE
 MAX. LOAD - 140 KG

REGULATION (EU) 2017/745 - CLASS 1 MD

FS34108
 FLADDERMUS | RESCUE STRETCHER

SPH236.01
 SPH FOR RESCUE HARNESS VERT 1

FS236.01
 VERT 1 - RESCUE HARNESS

EN361 + EN1497
 MAX. LOAD - 140 KG

REGULATION (EU) 2017/745 - CLASS 1 MD
 REGULATION (EU) 2016/425


(FS34109 = FS236.01 + SPH236.01)

FS34109
 BO | IMMOBILIZER VEST
 WITH INTEGRATED RESCUE HARNESS

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
FS34109 - BO | IMMOBILIZER VEST WITH INTEGRATED RESCUE HARNESS

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(EN)

INSTRUCTION MANUAL

WARNING: READ ALL INFORMATION CONTAINED IN THE TWO INSTRUCTIONS: GENERAL AND SPECIFIC.

SPECIFIC INSTRUCTIONS

Before using the rescue stretchers, you must read carefully and understand the safety information described on general instructions and the specific equipment instructions. The User Manual is an integral part of the device and therefore must be kept for the entire life of the device and must accompany it in any changes of use or ownership. ATTENTION!!! If you have any doubt about the product, if you need other language version of the instructions for use, declarations of conformity or any question about the product, please contact us: www.fallsafe-online.com. WARNING: The manufacturer and the seller decline any responsibility in case of incorrect use, improper application or modifications/repairs by persons not authorized by FALL SAFE®, invalidating the warranty and invalidating compliance with EU Regulation 2017/745. MODELS: FS34108 – FLADDERMUS I RESCUE STRETCHER and FS34109 – BO I IMMOBILIZER VEST WITH INTEGRATED RESCUE HARNESS (these models may be subject to implementation or change without notice). INTENDED USE AND CLINICAL BENEFITS: Rescue stretchers are devices intended to be used when the morphological characteristics of the intervention area require particularly compact equipment and adherence to the patient. FLADDERMUS and BO can also be lifted with the aid of hoisting systems fixed to the ground, using the integrated harness with the immobilizer vest. Patient intervention in the devices is not planned. PATIENTS: There are no particular indications related to the patient group. The product configuration is able to accommodate any subject as long as she/he is within the maximum capacity and within the limits of the size of the device. Expected patients are those who have injuries that prevent them from walking or are in an unconscious state, in a given rescue situation. CONTRAINDICATIONS, RISKS AND UNWANTED SIDE EFFECTS: No particular contraindications, residual risks or side effects are known with relation to the use of the device, as long as it is used in accordance with the users manual. USERS AND INSTALLERS: Please read section 1 – PHYSICAL CONDITIONING AND TRAINING from GENERAL INSTRUCTIONS. The intended users are rescue teams, carrying out operations related to the use of harnessing systems: 1 - Personnel trained for use of the device must also have training in managing lifting and handling suspended loads with people. 2 - Personnel who carry out interventions in situations classified as high risk or which are purely technical must be suitably trained and experienced in rescue. These devices are not intended for lay people. Do not allow untrained persons to help while using the product, as they may cause injury to themselves or others. The suitability of the users for use of this product can be attested by the training registration, in which trained persons, trainers, date and place are specified. This documentation must be kept for at least 10 years after the end of the product's life, and must be made available to the competent authorities and/or the manufacturer when requested. FALL SAFE® is always available for training. INSTALLATION: Installation is not required. STANDARDS: UE Regulation 2017/745 (UE Regulation on medical devices) As Distributor or End-User of the products manufactured and/or marketed by FALL SAFE®, users are strictly required to be familiar with the legal provisions in force in the country of destination of the goods, applicable to the devices to be supplied (including regulations relating to technical specifications and/or safety requirements) and, therefore, to understand the requirements necessary to ensure compliance of the products themselves with all legal requirements of the territory. With reference to EU Regulation 2017/745, please note that public or private operators who, when exercising their activity, detect an incident involving a medical product are required to notify the competent authority of the Member State, within the terms and in the manner established by one or more ministerial decrees, and notify the Manufacturer. MARKING AND TRACEABILITY: Each device is provided with a label, placed on the device itself, which contains the Manufacturer's identification data, product information, CE marking and serial number (SN). This must never be removed or covered, or else the warranty will be void as the device can no longer be traced. If the assigned SN cannot be traced, the device must be reconditioned, provided only under the responsibility of the manufacturer. EU Regulation 2017/754 requires manufacturers and distributors of medical devices to keep track of their location. If the device is in a location other than the address to which it was shipped or sold, or if it was donated, lost, stolen, exported or destroyed, permanently removed from use, or if the device was not delivered directly from FALL SAFE® please contact us or use FALL SAFE INSPECTOR®. WARNINGS: Use of the product for any purpose other than that described in the User Manual is prohibited. When using the devices, avoid contact with sharp or abrasive objects. The application of the device should not last longer than the time required for first aid operations and subsequent transport to the nearest rescue point. In case of exposed and/or injured skin, cover the surfaces in contact with the patient with a surgical sheet that respects bio-compatibility regulations to protect the patient's health. Always comply with the maximum capacity indicated in the User Manual. When determining the total weight load on the product, the operator should consider the weight of the patient, equipment and accessories. In addition, the operator should assess whether the overall size of the patient reduces the functionality of the product. Use only original or FALL SAFE® approved components/replacement parts and/or accessories to carry out any operation without causing alterations or modifications to the device. Before lifting, make sure that operators have a secure grip on the device. The device is a patient transport stretcher and cannot be used as a stationing device. For patient loading techniques for particularly heavy patients, for operations on steep terrain or in special and unusual circumstances, the presence of more operators is recommended in addition to the minimum foreseen operators. Always immobilize the patient. A failure to do so can cause serious damage to the patient. Do not move the stretcher if the weight is not properly distributed. Pay the utmost attention to any obstacles (water, ice, debris, etc.) present on the route, as they may cause the operator to lose balance and compromise proper functioning of the device. Dragging the stretcher on any type of surface leads to premature deterioration of the stretcher, which reduces its useful life and initial safety conditions. To preserve the life of the device, protect it as much as possible from UV rays and adverse weather conditions. Rescue stretchers are not approved for use with aircraft. BEFORE USE: Before each use, always check the conditions of the product, as specified in section 2 – BEFORE USE from "GENERAL INSTRUCTIONS". For first use, please check: if the packaging is intact and has protected the device during transportation; if all parts included in the packing list are present; general functionality of the device; product cleanliness; if there are no cuts, holes, lacerations or abrasions on the entire structure of the device and its parts (including detachable parts); correct fixing and holding of strap and ropes and the conditions of wear of the device and its standard parts. If the above conditions are met, the device may be considered ready for use. Non-compliant devices must be taken out of service. OPERATION: Before moving the patient, a medical evaluation of the patient's condition must be performed to stabilize the patient's clinical condition, verify possible surrounding hazardous situations, assess how to move the patient and whether it may be necessary to use alternative devices. Remove the device from the bag. Unlock all the straps and unroll the stretcher completely, while holding one end of it firmly on the ground with one foot. Then roll up the upper half of the stretcher in the opposite direction to that in which it was taken out of the bag, in order to eliminate the "memory" effect of the material. Carry out the same operation on the opposite side in order to make the stretcher perfectly flat. Place the stretcher on the ground. You can load the patient in two ways. Place the device next to the patient. Make sure that the side of the device on which the patient's head is to be placed is adjacent to the head itself. Position the transverse straps under the stretcher. Roll the patient to one side and slide the device as far as possible under the patient. Gently lay the patient on the device. Or fold the foot end of the stretcher to form a kind of ramp and place the belts under or to the side of the stretcher. Using

forearms to stabilize the head, grasp the patient in the underarm area, and then pull the patient onto the stretcher until the stretcher is in the proper position for anchoring. Make sure that you have correctly centred the patient in relation to the stretcher. Adjust all the harness straps around the patient's body and wrap her/him around the immobilizer vest. Raise the sides of the stretcher and fasten the straps, wrapping the stretcher around the patient. Do not cross the straps. Adjust the feet straps (highlighted in green) around the patient foot to prevent she/he from sliding. If only using BO, lift the patient by the D-ring with the A marking located at the chest (A). If only using FLADDERMUS, use the O-rings located at the sides of the stretcher; we recommend attaching FS806 to each side upper O-rings (O1) and to the anchor plate and attaching two anchor slings (FS800) to each side lower O-rings (O2), connecting FS910 to the slings and to the anchor plate, to facilitate changing the angle of the rescue stretcher; you can use the middle O-rings for extra anchor points. If using both, BO must be fixed to FLADDERMUS by the fix locks and by the rope loops (P2), then lift the patient by the rescue stretchers anchor points. The loops marked as P1 serve as a backup in the event of any anchor point becoming badly connected or failing. Transport by hand can only be performed if the patient is correctly positioned and all requirements specified in this user manual are met. Transporting the stretcher requires a minimum of 4 operators. A greater number of operators is necessary in the case of transport for long and/or difficult journeys. The handles are located on the sides of the rescue stretcher, attached to the straps (highlighted in orange). If it is not possible to transport the stretcher by hand, but ropes, winches, ladders, lifting straps or other suitable means are used, this is considered high risk rescue and must be carried out only by personnel adequately trained in this area and with the necessary experience on the field. MAINTENANCE, CLEANING AND INSPECTION: Please read "GENERAL INSTRUCTIONS. Checks must be made before and after use and regularly to verify: general functionality of the device, cleanliness of the device, holding and conditions of all straps and buckles, conditions of grommets, conditions of the rope knots and fulfillment of the requirements of the user manual. The frequency of inspections is determined by factors such as legal requirements, type of use, frequency of use, and environmental conditions during use and storage. All maintenance and overhaul activities must be recorded and documented with the relevant technical operation reports. This documentation must be kept for at least 10 years after the end of the product's life, and must be made available to the competent authorities and/or the Manufacturer when requested. Be sure to take all precautions to ensure that there is no risk of cross-infection or contamination of patients and operators. STORAGE: After carrying out cleaning and maintenance, lay out the rescue stretcher and place the straps facing inward; starting at the head end, roll up the stretcher as tightly as possible; continue rolling up the device, using your knees to prevent it from rolling out; fasten the straps to the metal buckles and place the device in its bag. For the immobilizer vest, fold down the sides and insert the device into its bag. Check section 8 – STORAGE/TRANSPORT from "GENERAL INSTRUCTIONS". WARRANTY AND DISPOSAL: Please read "GENERAL INSTRUCTIONS. ACCESSORIES AND SPARE PARTS: The accessories of these products consist of FS800 - ANCHOR SLING, FS806 - RESCUE Y-PAD and FS910 - EVACUATION KIT. We recommend the use of FS806 to avoid squeezing the victim when lifting her/him. FS910 is recommended for a faster and more versatile use of the rescue stretcher. This device has an adjustment mechanism that allows the rescue stretcher to be tilted to any desired angle. The connectors to be used can be selected considering the customer's needs and field of application. Spare parts can be supplied by the manufacturer. ATTENTION: ALSO READ THE MANUAL USER FOR FS236.01 and FS910.

ADDITIONAL INFORMATION

A11 - Acceptable Temperature; A12 - Storage; A13 - Annual Inspection; A14 - Cleaning; A15 - Drying; A16 - Hazards; A17 - Risk of death; A18 - Attention; A19 - Right; A110 - Wrong; A111 - Check.

MARKING/ LABELS

ML(A) - Serial Number; ML(B) - Reference number; ML(C) - Standard; ML(D) - Manufacture date; ML(E) - Brand name; ML(F) - Notified body which carries out product quality management; ML(G) - Instruction; ML(H) - QR code; ML(I) - CE Marking - complies with regulation (EU); ML(J) - Size; ML(K) - Maximum load; ML(L) - Fit icons.

EQUIPMENT RECORD

1-Product; 2-Reference Number; 3-Serial Number; 4-Manufacturing Date; 5-Purchase Date; 6-Date of first use; 7-Other relevant information; 8-Date; 9-Reason for entry; 10-Defects,Repairs, Etc; 11-Name & Signature; 12-Next periodic examination

NOMENCLATURE/FIELD OF APPLICATION

NFA1- Polyester Webbing - 45mm with Velcro straps®; NFA2- Polyester + Polyvinyl chloride - Textile; NFA3- EN1497 - Aluminium D-ring; NFA4- Polyester Webbing - 45mm; NFA5- IZY buckle - Aluminium; NFA6- EN361 - Aluminium D-ring; NFA7- Steel alloy "O" ring for webbing adjustment; NFA8- Plastic webbing keepers; NFA9- Polyester heavy duty seams; NFA10- Stainless Steel eyelet; NFA11- Static rope - Ø 9mm I Core: polyamide/ Cover: polyester; NFA12- Aluminium "O" ring for lift points; NFA13- Handle points - polyester webbing 45mm; NFA14- Stretcher body - High-density polyethylene; NFA15- Victim feet positioning - polyester webbing 45mm; NFA16- Polyester Webbing - 30mm; NFA17- Bottom rope loop - Static rope (11); NFA18- Steel alloy buckle - 3 bars; NFA19- Elastic webbing keepers; NFA20- Upper rope loop - Static rope (11).

P1/ P2 - See "SPECIFIC INSTRUCTIONS"

DONNING AND SETUP

FS34108

DS(1) - Buckle locations; DS(O) - Anchor points;

FS34109

DS(1) - Buckle locations; DS(A) - Anchor points; DS(P1) - Backup points; DS(P2) - FS34109 backup points and connection between FS34108 and FS34109

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